

Double Electric Breastpump

CALYPSO DOUBLE PLUS



Express milk in half the time

The Calypso Double Plus breastpump allows you to gently express milk from both breasts at the same time. This double expression helps increase and maintain the milk volume. The three different sizes of breast shell ensure that every mother can find the right size to suit her nipple. The soft Optiflow massage insert actively massages the breast and ensures efficient and gentle expression.



Key features at a glance

- Very quiet
- Simple, intuitive operation
- Extended guarantee
- The suction strength (vacuum) and frequency settings (cycle) can be adjusted individually and independently from each other
- Fulfills the highest standards of hygiene with its "Vacuum Seal" Technology
- Includes a cleaning brush
- Reduces the expression time by half, e.g. to 10 – 15 minutes
- Comes with three different breast shell sizes and an Optiflow massage insert





```
_V_("videojs11677", {"controls": true, "autoplay": false, "preload": false}).ready(function(){ var  
Playervideojs11677 = this; Playervideojs11677.volume(0.8); });
```

Related products: