

Low vacuum pump

Senator low vacuum



Need a reliable partner for performing demanding procedures involving low vacuum and reduced suction? Try the Ardo Senator low vacuum.

Features and Benefits:

- Sensitive – 15 litres per minute – the most delicate among suction pumps
- Fast – vacuum build-up to 100mbar
- Safe – incorporates “3guard” protective technology
- Whisper quiet – 37 dBA even at maximum vacuum
- Low maintenance – incorporates “LowFriction“
- Piston/Cylinder-technology
- Hygienic – easy clean surface and design
- Low vibration – completely still when in use

Specialist fields of use:

- Neurosurgery
- Pleural Drainage
- Various applications in secretion suction