

The Breastpump for professional use

CARUM




The new standard in Breastpumps for use in hospitals

With the Carum we are setting important standards and offer a professional quality hospital breastpump which meets the most demanding requirements of mothers and professionals.

[Get a quote now](#)

[Download brochure](#)

[Download case studies](#)



Important features at a glance

- Fully automatic expressing
- Vacuum and cycles are individually adjustable at any time, both in stimulation and expressing-mode
- Total safety and hygiene thanks to "Vacuum Seal" Technology
- "Sensitive Programme" in case of painful, inflamed nipples
- Especially gentle and comfortable transition between modes
- Double piston creates equal suction in both channels
- No additional consumables need replacing in or on the pump
- Ergonomic design makes the pump easy to clean
- Generous colour display with intuitive menu-driven management



```
_V_("videojs4541", {"controls": true, "autoplay": false, "preload": false}).ready(function(){ var  
Playervideojs4541 = this; Playervideojs4541.volume(0.8); });
```

Related products: